

PRODUCT INFORMATION SHEET

SPECIFICATIONS:

- 22"-28" Hydraulic Adjustable Seat Height
- 18" Stainless Steel Foot Ring
- 16" Diameter Round Seat
- 400 lb Weight Capacity

Overall Dimensions

Diameter of Seat	Adjustable Seat Height	Weight Capacity
16"	22"-28"	400 lbs

FEATURES:

- 3" Thick, 16" diameter round seat with high density foam
- Hands-free seat height range of 22" to 28"
- 18" Stainless Steel foot ring that is bolted
- Heavy-duty 3 1/2" thick hydraulic pump with a seat rotation brake that locks when foot pedal is lifted with heel of your shoe
- 2" Dual polyurethane casters
- Optional locking casters available (Item# LC)

FINISH SELECTIONS:

To view our **Vinyl Color Selection**, please visit:

http://www.swmedsource.com/Color%20Swatch%20pages/medical_vinyl_colors.html

CERTIFICATION(S):

As an option we offer CAL-133 protection, a flammability test procedure for seating furniture for use in Public Occupancies that has been adopted by many states. The test results show that the product passed both Criteria Q and Criteria B for the California TB-133 Burn Test.



COUNTRY OF ORIGIN:
United States of America



Provided by:

SW Med Source



Surgeon's Stool

Item #: SW-4080-R



Surgeon's Stool (with backrest)

Item #: SW-4080-R-B

WARRANTY INFORMATION:

All MC Stools have a 2 year warranty against defects in materials and workmanship.